

SPECIALTY & PE EQUIPMENT Optional Traveling Kits

MOVES MAT AT CEILING LEVEL



- Move your mat from wrestling area to adjoining storage area above the gym floor by installing a Traveling Kit. Mat travels at ceiling height, either end first or sideways, to the optimal position and then lowers to the floor.

FEATURES:

- Mat Lifter travel kit is driven by a ½ HP motor operating with the choice of 208-, 230- or 460-volts in 3-phase attached to a 60:1 ratio gear box with non-pneumatic wheel(s).
- Motors, drive assembly and all control mechanisms are encased in a fireproof cover.
- **End First Travel Kit** has two 1-ton machined trolleys with festooning power cables.
- **Sideways Travel Kit** has four half-ton machined trolleys with festooning power cables.

SUPPORTING DOCUMENTS

Access dimensions, specifications, instructions, technical drawings and other supporting documents located at:
draperinc.com/documents.aspx

For more information on this product visit:
draperinc.com/gymequipment/specialtype.aspx